

Performing by design

Sirius CTO

PTCA Balloon Catheter



Designed to cross the tightest lesions



ARTHESYS

Committed to excellence in developing Arterial Therapeutic Systems

Sirius CTO

PTCA Balloon Catheter

- > Unsurpassed crossability and pushability.
- > Lowest crossing profile.
- > Innovative shaft material for OTW version.


		Diameter (mm)			
		1.00	1.10	1.25	1.50
Length (mm)	8	□	□	□	□
	10	□	□	□	□
	15	□	□	□	□
	20	□	□	□	□
	25	□	□	□	□

		Pressure (bars)									
		2	4	6*	8	10	12	14	16	18**	20
Balloon diameter (mm)	ø1.00	0.96	0.98	1.00	1.02	1.04	1.06	1.08	1.10	1.12	1.15
	ø1.10	1.06	1.08	1.10	1.12	1.14	1.17	1.19	1.21	1.23	1.27
	ø1.25	1.20	1.23	1.25	1.28	1.30	1.32	1.35	1.38	1.40	1.44
	ø1.50	1.44	1.47	1.50	1.53	1.56	1.59	1.62	1.65	1.68	1.73


□ Available

* Nominal pressure ** Rated Burst Pressure

Sirius CTO RX



Sirius CTO OTW



- > Guidewire compatibility: 0.014" maximum recommended
- > Guiding catheter compatibility: 5F (0.058" ID)



Take a look at our most up-to-date developments and latest news
info@arthsys.com
<http://www.arthsys.com>

20, rue Traversière
 92230 GENNEVILLIERS - FRANCE
 Tel: + 33 (0)1 4111 8777
 Fax: + 33 (0)1 4111 8770